



Partnering For Success

QS/1<sup>®</sup>

| October 2009

INSIGHT

# QS/1 &



# PickPoint

PickPoint's WCS interfaces to the QS/1 Pharmacy Management System for fast, accurate, safe, compliant and efficient Will Call management. Utilizing leapfrog technology, your pharmacy will build confidence, professionalism and sales with your patients and staff.

With the unique interface to QS/1's powerful pharmacy management systems, all necessary information is transferred from QS/1 systems to WCS. WCS provides fail-safe backup; all data is secure and does not leave the control of the pharmacy.

The prescriptions are placed in a clear bag affixed to an intelligent hanger and can be placed randomly on an organized wall rack. The rack is scalable in size and is produced by the nationally known and respected UniWeb fixture company. The hanger bars on the wall are moveable and serve to provide low voltage charging capability to the hanger bags which have internal lithium batteries for quick charge and long life.

It is important to note that each hanger/bag, bar and wall rack is unique, thus preventing any single point of failure that might otherwise interrupt the system at that most important point when a customer is at the counter. Random placement of bags is not only fast but efficient from a space and appearance standpoint. Most pharmacies reduce the amount of space previously required by their will call storage bins.

When a patient comes to the pickup counter, after identification in any one of a number of optional ways, the bag, or bags in the case of multiple prescriptions being submitted or completed for the patient, will light up on the patented hanger/bag. Each pickup workstation is user-assigned a light color that illuminates for pulling. If the wrong bag is selected, it will not clear for release.

Both Return to Stock and Interactive Voice Response functions are facilitated by WCS. It has a parameter-setting capability of lighting or electronically messaging user-defined cycle or deadline dates are approaching or have been reached. RTS becomes a daily function by lighting those prescriptions which are to be pulled. The pharmacy's IVR is queued by WCS to place a patient call when the Rx is ready for pickup.

Go to [www.PickPoint.com](http://www.PickPoint.com) for more information.

## Your Rx is Ready for Pickup

by Tom George  
Sr. Vice President Pharmacy Services  
PickPoint Corporation

